



CYCLONE® KIUBE

—  
Max potential,  
max capacity.

RADIO  
PHARMA  
SOLUTIONS

# 35 years of experience concentrated in one innovation.



## — 30 years of innovation and expertise concentrated in the CYCLONE®KIUBE

### The first world reference

In 1986, IBA launched the first compact cyclotron for industrial production of radioisotopes. This revolutionary new design increased the production efficiency by a factor 15 and created a world reference.

30 years later, IBA once again makes history with the launch of its new 18 MeV cyclotron. CYCLONE®KIUBE is more compact and the most powerful mid energy PET cyclotron. It has been developed with three keywords in mind: **Reliability, High Performance and Flexibility**. These keywords were defined in collaboration with our Users.

### Expertise benefiting to a large installed base

With IBA as partner, you have the assurance to remain at the forefront of innovation and to benefit of every novel feature for tracer production.

More than 300 cyclotrons installed over 5 continents demonstrate this expertise and customer satisfaction.

### Support over the lifetime of your system

IBA provides continuous support over the lifetime of your project, expanding your skills and applications, boosting your uptime and maximizing the return on your investment with all the latest innovations, services and training.

In addition, the system is built so that future upgrades can be easily implemented to keep your cyclotron state-of-the-art.


# 150

IBA PATENTS  
ON PARTICLE  
ACCELERATORS

# 35

YEARS OF  
EXPERTISE

## — A story of Cyclotron Innovations



# Designed for Ever.

— CYCLONE®KIUBE will evolve with you, for you

## Upgradable like no other

CYCLONE®KIUBE is upgradable like no other cyclotron, so you can increase your production capacity; from 100 µA to 150µA, 180 µA or even up to 300µA on target.

With a lower initial investment, you can start with the CYCLONE®KIUBE 100 while maintaining the ability to increase your PET Center's capacity over time with an on-site upgrade.

## Infinite flexibility


With eight independent exit ports, CYCLONE®KIUBE is the most flexible system, producing the widest range of PET radioisotopes.

A full range of Nirta® targets are available giving you access to <sup>18</sup>F, <sup>13</sup>N, <sup>15</sup>O, <sup>11</sup>C (CO<sub>2</sub> & CH<sub>4</sub>), <sup>18</sup>F<sub>2</sub>, <sup>68</sup>Ga, ...

In addition, IBA provides Nirta® Solid target technology for the production of radioisotopes such as <sup>64</sup>Cu, <sup>89</sup>Zr, <sup>124</sup>I, ...

## Custom Energy.

The IBA-patented custom energy CYCLONE®KIUBE allows users to irradiate the target directly at a fixed energy level of 13, 14 or 15 MeV, without using a degrader.



This is our typical production schedule thanks to the IBA 8-exit ports cyclotron:

Time	Radionuclide	Radio pharmaceutical
4.30-6.00	<sup>18</sup> F	FDG
6.30-8.00	<sup>18</sup> F	F-other
8.15-8.45	<sup>11</sup> C	CHOLINE
10.00-12.00	<sup>18</sup> F	MK9470
12.30-13.00	<sup>11</sup> C	Raclopride
13.30-13.40	<sup>13</sup> N	NH3
14.00-16.00	<sup>15</sup> O	10*H <sub>2</sub> O
16.30-17.00	<sup>11</sup> C	PDE10A



**Prof. Dr Guy Bormans**  
Head of radiopharmacy  
KATHOLIEK UNIVERSITEIT LEUVEN, BELGIUM



### CYCLONE® KIUBE WITH SELF-SHIELDING OPTION\*



**10Ci**  
370 GBq  
100 µA

**16Ci**  
592 GBq  
150 µA

UPGRADABLE

### CYCLONE® KIUBE VAULTED\*



**10Ci**  
370 GBq  
100 µA

**16Ci**  
592 GBq  
150 µA

**20Ci**  
740 GBq  
180 µA


**30Ci**  
1110 GBq  
300 µA

UPGRADABLE

\*Production in Curies of <sup>18</sup>F- in 2 hours run



# Designed to Deliver.



## After 10 years of excellent experiences with Cyclone® 18/9,

*we have now added the new CYCLONE®KIUBE. Outstanding! All well thought engineering details made the operation unprecedentedly easy, reliable and maintenance fast and safer, lowering dose exposure.*



**Francisco Alves**  
Chief physicist & head of Cyclotron  
ICNAS-UNIV. COIMBRA - PORTUGAL



# 300 $\mu\text{A}$

MAX CURRENT OF THE CYCLONE®KIUBE

## — CYCLONE®KIUBE delivers outstanding performances

### The most cost-effective

It has been demonstrated that 18 MeV is the optimal proton energy for the highest production yield of most of PET radioisotopes. Your production capacity is optimized so as your revenue; your center footprint is kept to a reasonable size so as your investment. This means that 18 MeV is the most cost-effective solution for your PET isotopes production.

The growing demand for radioisotopes means a greater need for efficiency. CYCLONE®KIUBE saves enriched water and has the lowest power consumption per Curie produced, yet is the most powerful of the market.

### The most performant cyclotron

CYCLONE®KIUBE offers unmatched production capacity for an internal source PET cyclotron. A 2-hour dual beam run could generate up to 30 Ci of  $^{18}\text{F}$ -; hence, an incredible batch of FDG synthesized on your Synthera®+ modules.

### Maximum reliability


Its unique twin-proton sources give maximum reliability, as switching to the second source during operation is simple, fast and completely automated, translating into more than 99% source uptime. Moreover, the production is optimized in real time thanks to the automatic tuning of the ion source, stripper and magnetic field.

#### HIGHEST PRODUCTION CAPACITY

$\mu\text{A}$	$^{18}\text{F}$ - capacity/run	FDG doses/run[*]
CYCLONE®KIUBE 100	2 hours : 10 Ci / 370 GBq	60 - 100
CYCLONE®KIUBE 150	2 hours : 16 Ci / 592 GBq	96 - 160
CYCLONE®KIUBE 180	2 hours : 20 Ci / 740 GBq	120 - 200
CYCLONE®KIUBE 300	2 hours : 30 Ci / 1110 GBq	180 - 300

[\*] using Synthera® FDG chemistry (12h shelf life - stability).  
May vary according to local conditions, transport time and patients schedule.

# Designed for you.



**<100 M<sup>2</sup>**  
1,070 FT<sup>2</sup> PET CENTER AREA

## Compact design

Its compact design can fit into any existing cyclotron vault.

Cyclone®KIUBE is also available with a self-shielding option.

## IntegraLab® ONE

IBA experts designed the **world's most optimized** [c]GMP PET center with a footprint of less than 100m<sup>2</sup> (1,070 square feet) including a 18 MeV self-shielded cyclotron, while allowing a very high production capacity of <sup>18</sup>F for FDG and <sup>68</sup>Ga.

IntegraLab®ONE is a real ready-to-run integrated Radiopharmacy, designed to ensure a fast, smooth and risk-free set up of your [c] GMP radiopharmacy.

## IntegraLab® PLUS

IntegraLab® PLUS is a **fully integrated** solution, combining the equipment and services required to establish a Radiopharmaceutical [c]GMP Production Center tailored to your precise needs.

Our approach is **truly comprehensive**, taking your project from building design to production, with full regulatory compliance, and with the selection, integration, installation and qualification of the equipment needed to meet your radioisotope production goals.

## — A self-monitoring cyclotron that maximizes your uptime

### Unmatched reliability

CYCLONE®KIUBE has an unmatched reliability. Its compact design includes redundancy for the main systems (ion sources, strippers, targets, vacuum pumps,...).



The software self-tests prior to production and automatically self-audits after maintenance, leading to maximized uptime.

The choice of low activation material combined with the optimization of the maintenance tasks ensure the lowest radiation exposure for your staff.

### User-friendly Software

The software is very user-friendly, meaning it is easier to train and rotate the staff that operates the system. It provides an easy and guided control of the system with a simple production mode for an automatic set up to maximise your <sup>18</sup>F production and the manual mode gives you full flexibility.

The 24/7 helpdesk provides you with peace of mind, with IBA experts analyzing cyclotron data in order to diagnose and solve 95% of problems remotely.





# Nirta® target technology

Complete target range for your extended needs

## — <sup>18</sup>F- conical liquid target

### High power with low enriched water consumption

- Efficient design to achieve 4Ci/ml in 2 hour-run
- Lowest enriched water consumption on the market
- Wide range of volumes available: your production capacity will always match your needs
- Quick and low dose maintenance: less pieces and o-rings
- Easier disassembling, assembling and manipulation.



## <sup>18</sup>F

	Conical 5	Conical 8	Conical 12	Conical 16
Chemical form		[ <sup>18</sup> F-]		
Target reaction		<sup>18</sup> O(p,n) <sup>18</sup> F		
Target material		H <sub>2</sub> <sup>18</sup> O		
Window material		Havar®		
Insert material		Niobium		
Enriched water volume [ml]	1.8 - 2.0	2.3 - 2.5	2.7 - 3.0	3.5 - 4.3
Irradiation time [min]	120	120	120	120
Recovered activity EOB [Ci] / [GBq]	5 / 185	8 / 296	12 / 444	16 / 592

Thanks to the combination of the CYCLONE®KIUBE and the Nirta® <sup>18</sup>F- targets,

we have a cost-effective solution for large production of <sup>18</sup>F; high activity in a short run with low enriched water consumption. We are happy with the reliability and performance and we are proud to produce 60 Ci <sup>18</sup>F a day to deliver to hospitals for patient diagnosis.



**Volkan Çam**  
Technical Manager  
NUKLEON, TURKEY

## — Liquid targets

### Unique solution for <sup>68</sup>Ga production

The unique and patented [US10600528B2 and EP3101660B1] solution for production of <sup>68</sup>Ga-radiopharmaceuticals for human use has many advantages :

- Seamless automation of the whole process with the Synthera® platform
- No risk of long-lived <sup>68</sup>Ge [271 days] in the final product
- Production every hour depending on patient schedule
- Simple & cost effective GMP production process for in-house use or for distribution production similar to the <sup>18</sup>F process



	<sup>68</sup> Ga	<sup>13</sup> N	<sup>64</sup> Cu <sup>(1)</sup>
Chemical form	[ <sup>68</sup> Ga]GaCl <sub>3</sub>	[ <sup>13</sup> N]NH <sub>3</sub>	[ <sup>64</sup> Cu]CuCl <sub>2</sub>
Target reaction	<sup>68</sup> Zn(p,n) <sup>68</sup> Ga	<sup>16</sup> O(p,α) <sup>13</sup> N	<sup>68</sup> Ni(p,n) <sup>68</sup> Cu
Target material	Enriched <sup>68</sup> Zn salt	H <sub>2</sub> O (natural) + 5mMol Ethanol	enriched <sup>64</sup> Ni solution
Window material	Niobium	Havar®	Niobium
Insert material	Niobium	Niobium	Niobium
Energy degrador	na	Graphite	na
Beam energy	18 MeV	16 MeV	18 MeV
Target yield [mCi/μA sat]	~(2)	38	~(2)
irradiation time [min]	60	20	300
Insert volume [ml]	3.7	3.7	5
Recovered activity EOB [mCi] / [GBq]	200 / 7.4	1000 / 37	50/1.85

<sup>(1)</sup> not commercially available

<sup>(2)</sup> depending on target solution concentration and other parameters

## — Gas targets

	<sup>15</sup> O	<sup>11</sup> C – CO <sub>2</sub>	<sup>11</sup> C – CH <sub>4</sub>
Chemical form	[ <sup>15</sup> O]O <sub>2</sub>	[ <sup>11</sup> C]CO <sub>2</sub>	[ <sup>11</sup> C]CH <sub>4</sub>
Target reaction	<sup>15</sup> N(p,n) <sup>15</sup> O	<sup>14</sup> N(p, α) <sup>11</sup> C	<sup>14</sup> N(p, α) <sup>11</sup> C
Target material	<sup>15</sup> N <sub>2</sub> + 0.5-1% O <sub>2</sub>	N <sub>2</sub> + 0.5-1% O <sub>2</sub>	N <sub>2</sub> + 5% H <sub>2</sub>
Window material	Aluminum	Aluminum	Aluminum
Insert material	Aluminum	Aluminum	Aluminum
Energy degrador	Graphite	na	na
Beam energy	13 MeV	18 MeV	18 MeV
Insert volume [ml]	4 – 5/min (cont. flow)	55	40
Target yield[mCi/μA sat]	30	150	60
irradiation time [min]	online	30	30

## — Solid targets

### Solid target technology for the production of novel radioisotopes

- Routine production and research programs
- Full packages available including dissolution and purification modules
- Possibility to have a pneumatic transport system between the cyclotron vault and the radiochemistry laboratory.
- The Nirta® Plus solid target can be installed inside the self-shielding



Nirta® solid target



Nirta® Plus solid target



Nirta® High Power<sup>(1,2)</sup> solid target

Target current	20-40 μA	40-75μA	up to 300μA <sup>3</sup>
Target angle	90°	35°	customizable capsule (from 10°-90°)
Automated transfer	1 way in option	1 way in option	2 ways

<sup>1</sup> EP18187472.8.

<sup>2</sup> Vectio beam line required

<sup>3</sup> Energy degrader on capsule will reduce beam current capacity to 150 μA

### TYPICAL PRODUCTION WITH SOLID TARGETS – PET/SPECT ISOTOPES

	<sup>64</sup> Cu <sup>(1)</sup>	<sup>89</sup> Zr <sup>(1)</sup>	<sup>124</sup> I <sup>(1)</sup>	<sup>123</sup> I <sup>(1)</sup>	<sup>66</sup> Ga	<sup>67</sup> Ga	<sup>111</sup> In	<sup>86</sup> Y
Half-life	12.7h	78h	4.18d	13.27h	9.49h	78h	67h	14.7h
Target reaction	<sup>64</sup> Ni(p,n) <sup>64</sup> Cu	<sup>89</sup> Y(p,n) <sup>89</sup> Zr	<sup>124</sup> Te(p,n) <sup>124</sup> I	<sup>123</sup> Te(p,n) <sup>123</sup> I	<sup>66</sup> Zn(p,n) <sup>66</sup> Ga	<sup>67</sup> Zn(p,n) <sup>67</sup> Ga		<sup>86</sup> Sr(p,n) <sup>86</sup> Y
Target Material	<sup>64</sup> Ni	<sup>89</sup> Y	<sup>124</sup> Te	<sup>123</sup> Te	<sup>66</sup> Zn	<sup>67</sup> Zn	<sup>111</sup> Cd	<sup>86</sup> Sr
Energy [on target] [MeV]	14	15	15	16	16	16	16	18
Yield [mCi/μAh] <sup>(2)</sup>	1.34-8.5 <sup>(3)</sup>	0.34	0.45	2.34	-	-	-	-

<sup>(1)</sup> commercial package available

<sup>(2)</sup> The production rate is given by yield x target current x irradiation time.

The maximum target current depends on the irradiation station (Nirta® Solid, Nirta® Solid Plus, Nirta® High power solid)

<sup>(3)</sup> depending on target material thickness. 0.073 mCi/μAh/mg of Ni-64

Enlarge your cyclotron possibilities.

The Vectio® beam line extension can be installed for high power solid target work and for research with proton beam.

# Technical features.



	Cyclone® KIUBE 100	Cyclone® KIUBE 150	Cyclone® KIUBE 180	Cyclone® KIUBE 300
<b>HIGH CAPACITY PROTON BEAM</b>				
Energy	18 MeV			
Beam current on target	100 µA	150 µA	180 µA	300 µA
<b>TARGET FLEXIBILITY</b>				
Number of target ports	8			
Simultaneous extracted beams	2			
Vectio® Beam Transport Line	short, medium, long			
<b>LOW POWER CONSUMPTION</b>				
Stand-by mode	< 3 kW			
Beam-on mode	< 45 kW	< 55 kW	< 60 kW	< 65 kW
<b>COMPACT DESIGN</b>				
Cyclotron weight	18 Tons			
Cyclotron overall dimensions [m]	1.9 x 1.9 x 1.8 (l x w x h)			
Internal room dimensions [m]	3.8 x 3.8 x 2.5 - Recommended : 4 x 4 x 2.5 (l x w x h)			
Internal room dimensions with self-shielding [m]	6 x 7 x 3 (l x w x h)		-	-
Self-shielding overall dimensions/piece [m]	4.645 x 2.770 x 2.575 (l x w x h)		-	-

## ABOUT IBA (Ion Beam Applications S.A)

IBA is a cancer diagnostics and treatment company and the worldwide technology leader in the field of proton therapy. The company's expertise lies in the development of next-generation proton therapy technologies and radiopharmaceuticals that provide oncology care providers with premium quality services and equipment, including IBA's leading fully-integrated IntegraLab® system.

## ABOUT IBA RADIOPHARMA SOLUTIONS

Based on long-standing expertise, IBA RadioPharma Solutions supports hospitals and radiopharmaceutical distribution centers with their in-house radioisotope production by providing them with global solutions, from project design to the operation of the facility. In addition to high-quality technology production equipment, IBA has developed in-depth experience in setting up GMP radiopharmaceuticals production centers.



### CONTACT US:

Chemin du Cyclotron, 3 | 1348 Louvain-la-Neuve, Belgium  
Tel.: + 32 10 201 275 | Fax: + 32 10 47 59 58 |  
RadioPharmaSolutions@iba-group.com